

Product components

Components	Component number	Size		Size	
		20	RXN	100	RXN
10X ABScript II Enzyme Mix	RM21451	40	μL	200	μL
2X ABScript II Reaction Mix	RM21450	200	μL	1	mL
Oligo d(T) ₂₃ VN * (50 μM) **	RM20115	40	μL	200	μL
Random Primer Mix (60 μM) **	RM20116	40	μL	200	μL
dNTPs (10 mM each)	RM20120	30	μL	100	μL
Nuclease-free H ₂ O	RM20214	1.25	mL	1.25	mL

* V = A, G or C; N = A, G, C or T.

** Contains 1 mM dNTP.

Product Description

ABScript II cDNA First Strand Synthesis Kit features two optimized mixes, ABScript II Enzyme Mix and ABScript II Reaction Mix. The enzyme mix combines ABScript II Reverse Transcriptase and RNase Inhibitor, and the reaction mix contains an optimized buffer. ABScript II Reverse Transcriptase is a recombinant M-MuLV reverse transcriptase with reduced RNase H activity and increased thermostability. It can be used to synthesize first strand cDNA at higher temperatures than the wild type M-MuLV. The enzyme is active up to 48°C, which provides higher specificity and higher yield of cDNA.

The kit also provides two optimized primers for reverse transcription and nuclease-free water. An anchored Oligo-dT primer [d(T)₂₃VN] forces the primer to anneal to the beginning of the polyA tail. The optimized Random Primer Mix provides random and consistent priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs. The first strand cDNA product generated is up to 10 kb.

Storage

Store at -20°C

First Strand cDNA Synthesis Reaction

1. Denaturation of RNA and primer at 65–70°C for 5 minutes can remove secondary structures that may impede long cDNA synthesis. However, this step can be omitted in many cases (unpublished results).
2. We recommend incubation at 42°C for one hour for maximum yield and length. However, many targets can be detected after a much shorter incubation time. For example, 10 minutes incubation can be used for up to 5 kb cDNA synthesis.

Choice of Primers for Reverse Transcription

1. Oligo-dT primer is preferred for most applications because it ensures that all cDNA copies terminate at the 3' end of the mRNA and produces the longest contiguous cDNA. An anchored oligo-dT primer [d(T)₂₃VN] forces the primer to anneal to the start of the polyA tail, thereby preventing priming at internal sites in the polyA tail⁽¹⁾.
2. The Random Primer Mix is an optimized mix of hexamer and dNTP. It provides random priming sites covering the entire RNA templates including both mRNAs and non-polyadenylated RNAs (such as ribosomal RNAs). The Random Primer Mix yields shorter cDNAs on average and can be used for the detection of multiple short RT-PCR products. Random Primer Mix offers good performance in a wide range of RNA templates.
3. When a gene-specific primer is used in a cDNA synthesis reaction, the cDNA product can be used only for amplification of that transcript. This priming method gives good results when the amount of RNA is limiting (below 10 ng) and only one particular cDNA is desired.

4. Recommended primer concentration:

Primer	Oligo d(T) ₂₃ VN	Random Primer Mix	Specific Primer
Final Conc.	5 µM	6 µM	0.1-1 µM

First Strand cDNA Synthesis Protocols

Thaw kit components on ice and mix by inverting several times.

Easy Protocol

- Mix the following components and incubate at 42°C for 1 hour. If Random Primer Mix is used, an incubation step at 25°C for 5 minutes is recommended before the 42°C incubation.

Components	1X RXN
Template RNA	up to 1 µg
Oligo d(T) ₂₃ VN (50 µM) *	2 µL
10 mM dNTPs	1 µL
2X ABScript II Reaction Mix	10 µL
10X ABScript II Enzyme Mix	2 µL
Nuclease-free H ₂ O	to 20 µL

***Note:** According to different requirements, primers can be added according to the recommended primer concentrations: Oligo d(T)₂₃VN (5 µM), random primers (6 µM), and specific primers (0.1-1 µM).

- Inactivate the enzyme at 80°C for 5 minutes. For downstream PCR application, the volume of cDNA product should not exceed 1/10 of the PCR reaction volume.

Standard Protocol

Note: If denaturation of template RNA is desired, use the following protocol. Standard protocol helps to denature the secondary structure of complex RNA templates, improve the efficiency of reverse transcription, and increase the yield of cDNA products.

- Mix RNA sample and primer d(T)₂₃VN in a sterile RNase-free microfuge tube.

Components	1X RXN
Total RNA	1-5 µL (up to 1 µg)
Oligo d(T) ₂₃ VN (50 µM) *	2 µL
10 mM dNTPs	1 µL
Nuclease-free H ₂ O	to 8 µL

***Note:** According to different requirements, primers can be added according to the recommended primer concentrations: Oligo d(T)₂₃VN (5 µM), random primers (6 µM), and specific primers (0.1-1 µM).

- Denature sample RNA/d(T)₂₃VN for 5 minutes at 65°C. Spin briefly and put promptly on ice.
- Add the following components.

Components	1X RXN
2X ABScript II Reaction Mix	10 µL
10X ABScript II Enzyme Mix	2 µL

- Incubate the 20 µL cDNA synthesis reaction at 42°C for one hour. If Random Primer Mix is used, an incubation step at 25°C for 5 minutes is recommended before the 42°C incubation.
- Inactivate the enzyme at 80°C for 5 minutes. The cDNA product should be stored at -20°C. In general, the volume of cDNA product should not exceed 1/10 of the PCR reaction volume.

No-RT Negative Control Reaction

Mix the following components and incubate at 42°C for 1 hour.

Components	1X RXN
Template RNA	up to 1 µg
Oligo d(T) ₂₃ VN (50 µM)	2 µL
2X ABScript II Reaction Mix	10 µL
10 mM dNTPs	1 µL
Nuclease-free H ₂ O	to 20 µL

General Information for Successful cDNA Synthesis

1. Intact RNA of high purity is essential for sensitive RT-PCR detection.
2. Total RNA or mRNA can be used in the reverse transcription reaction. Total RNA is generally sufficient for most RT-PCR analyses.

Trouble-shooting Guide

Low Yield of cDNA

1. Check the integrity of the RNA by denaturing agarose gel electrophoresis ⁽²⁾.
2. RNA should have a minimum A260/A280 ratio of 1.7 or higher. Ethanol precipitation followed by a 70% ethanol wash can remove most contaminants such as EDTA and guanidinium. Precipitation with lithium chloride can remove polysaccharides ⁽²⁾.
3. Phenol/chloroform extraction and ethanol extraction can remove contaminant proteins such as proteases ⁽²⁾.
4. Some target RNA may contain strong pauses for RT; Use random priming instead of d(T)₂₃VN.
5. Use sufficient amount of RNA.

References

1. Liao, J. and Gong, Z. (1997) Biotechniques 23, 368-370.
2. Sambrook, J. and Russel, D.W. (2001). Molecular Cloning: A Laboratory Manual, (3rd ed.), (pp. 8.46-8.53 and 11.37-11.42). Cold Spring Harbor: Cold Spring Harbor Laboratory Press.