

Product Components

Components	Component number	Size-1	Size-2
DL1000 DNA Marker	RM30291	500 μ L	5 $\times$ 500 μ L

Product Description

The DL1000 DNA Marker is a ready-to-use DNA Marker, which is pre-mixed with loading dye for direct gelloading. This product consists of 8 individual DNA fragments: 1000 bp, 2000 bp, 3000 bp, 4000 bp, 5000 bp, 6000 bp, 8000 bp and 10000 bp, which are obtained by enzyme digestion of recombinant plasmids.

Storage

Store at -20°C for long-term storage. Store at 4°C after thawing for short-term use and avoid repeated freezing and thawing.

Notes

1. Please thaw and mix thoroughly before use.
2. Please change the electrophoresis buffer in time, and use freshly prepared agarose gel to avoid affecting the electrophoresis results.

Experiment Process

1. As a ready-to-use product, load 5 μ L on gel for each sample is recommended, and if the well is wide, loading volume can be increased appropriately for electrophoresis.
2. The recommended condition is 1 $\times$ TAE buffer or 0.5 $\times$ TBE buffer, 1% agarose gel, and 4-10 V/cm between positive and negative electrodes.
3. Use UV transilluminator to detect DNA bands stained with nucleic acids stain dye.
4. The electrophoresis results of the DL1000 DNA Marker are shown below.

