

PPP3R1 Rabbit pAb

Catalog NO.: A0954

Basic Information

Observed MW

15kDa

Calculated MW

19kDa

Category

Polyclonal Antibodies

Applications

WB, IF/ICC, ELISA

Cross-Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Background

Enables cyclosporin A binding activity; phosphatase binding activity; and protein domain specific binding activity. Involved in calcineurin-NFAT signaling cascade and positive regulation of transcription by RNA polymerase II. Part of calcineurin complex. Implicated in Alzheimer's disease and dilated cardiomyopathy.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

ELISA

Recommended starting concentration is 1 µg/mL. Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
5534	P63098

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CNB; CNB1; CALNB1; PPP3R1

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates using PPP3R1 Rabbit pAb (A0954) at 1:1000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.



Immunofluorescence analysis of NIH/3T3 cells using

PPP3R1 Rabbit pAb (A0954) at dilution of 1:100.

Secondary antibody: Cy3-conjugated Goat anti-Rabbit

IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear

staining.

