

CGN Rabbit pAb

Catalog NO.: A15489 **1 Publications**

Basic Information

Observed MW

160kDa

Calculated MW

137kDa

Category

Polyclonal Antibodies

Applications

WB, IHC-P, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Enables cadherin binding activity. Predicted to act upstream of or within bicellular tight junction assembly; epithelial cell morphogenesis; and microtubule cytoskeleton organization. Located in bicellular tight junction and plasma membrane.

Recommended Dilutions

WB 1:200 - 1:2000

IHC-P 1:50 - 1:200

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

57530

Swiss Prot

Q9P2M7

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CGN; cingulin

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using CGN Rabbit pAb (A15489) at 1:1000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma using CGN Rabbit pAb (A15489) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.

