

GPRIN1 Rabbit pAb

Catalog NO.: A15555

Basic Information

Observed MW

115kDa

Calculated MW

102kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Background

Predicted to enable phosphoprotein binding activity. Predicted to be involved in neuron projection development. Predicted to be located in growth cone. Predicted to be active in plasma membrane.

Recommended Dilutions

WB 1:200 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

114787

Swiss Prot

Q7Z2K8

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

GRIN1; GPRIN1


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at -20°C. Avoid freeze / thaw cycles.		
Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.		

Validation Data



Western blot analysis of various lysates using GPRIN1 Rabbit pAb (A15555) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.

