

ANP32E Rabbit pAb

Catalog NO.: A17220 **2 Publications**

Basic Information

Observed MW

31 kDa

Calculated MW

31 kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables histone binding activity. Involved in histone exchange. Located in nucleus. Part of Swr1 complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

81611

Swiss Prot

Q9BTT0

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

LANPL; LANP-L; ANP32E

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using ANP32E Rabbit pAb (A17220) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 10s.

