

TMF1 Rabbit pAb

Catalog NO.: A17542

Basic Information

Observed MW

160kDa

Calculated MW

123kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Conjugation

Unconjugated

Background

Enables androgen receptor binding activity and nuclear receptor coactivator activity. Involved in androgen receptor signaling pathway and positive regulation of transcription by RNA polymerase II. Located in Golgi apparatus.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

7110

Swiss Prot

P82094

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

TMF; ARA160; TMF1


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at -20°C. Avoid freeze / thaw cycles.		
Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.		

Validation Data



Western blot analysis of various lysates using TMF1 Rabbit pAb (A17542) at 1:1000 dilution.
Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 180s.

