

SLTM Rabbit pAb

Catalog NO.: A18493

Basic Information

Observed MW

117kDa

Calculated MW

117kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

Enables RNA binding activity. Predicted to be involved in regulation of mRNA processing and regulation of transcription by RNA polymerase II. Located in nuclear body.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

79811

Swiss Prot

Q9NWH9

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

Met; SLTM


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at -20°C. Avoid freeze / thaw cycles.		
Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.		

Validation Data



Western blot analysis of lysates from Mouse brain, using SLTM Rabbit pAb (A18493) at 1:1000 dilution.
 Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
 Lysates/proteins: 25µg per lane.
 Blocking buffer: 3% nonfat dry milk in TBST.
 Detection: ECL Basic Kit (RM00020).
 Exposure time: 30s.

