

PLSCR3 Rabbit mAb

Catalog NO.: A20915

Recombinant

1 Publications

Basic Information

Observed MW

32kDa

Calculated MW

32kDa

Category

Monoclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Clone number

ARC2887

Background

Enables several functions, including calcium-dependent protein binding activity; metal ion binding activity; and phospholipid scramblase activity. Involved in several processes, including cardiolipin biosynthetic process; regulation of apoptotic process; and regulation of release of cytochrome c from mitochondria. Located in cytosol; mitochondrion; and plasma membrane.

Recommended Dilutions

WB

1:500 - 1:1000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
57048	Q9NRY6

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

PLSCR3

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



