

DLGAP5 Rabbit pAb

Catalog NO.: A2197 **2 Publications**

Basic Information

Observed MW

120kDa

Calculated MW

95kDa

Category

Polyclonal Antibodies

Applications

WB, IF/ICC, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Background

Predicted to enable microtubule binding activity. Predicted to be involved in several processes, including centrosome localization; kinetochore assembly; and mitotic spindle organization. Located in several cellular components, including centrosome; cytosol; and mitochondrion.

Recommended Dilutions

WB 1:500 - 1:2000

IF/ICC 1:50 - 1:200

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
9787	Q15398

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

DLG7; HURP; DLGAP5

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of lysates from K562 cells, using DLGAP5 Rabbit pAb (A2197) at 1:1000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.



Immunofluorescence analysis of U2OS cells using DLGAP5 Rabbit pAb (A2197). Secondary antibody: Cy3-conjugated Goat anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.

