

ASUN/INTS13 Rabbit mAb

Catalog NO.: A22792 **Recombinant**

Basic Information

Observed MW

80kDa

Calculated MW

80kDa

Category

Monoclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human

Conjugation

Unconjugated

Clone number

ARC55758

Background

Involved in regulation of mitotic cell cycle. Acts upstream of or within centrosome localization; mitotic spindle organization; and protein localization to nuclear envelope. Located in cytoplasm and nuclear body.

Recommended Dilutions

WB 1:10000 - 1:80000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID	Swiss Prot
55726	Q9NVM9

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

ASUN; GCT1; NET48; Mat89Bb; SPATA30; C12orf11; ASUN/INTS13

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates, using ASUN/INTS13 Rabbit mAb (A22792) at 1:70000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 60s.

