

MTA3 Rabbit pAb

Catalog NO.: A2328 **2 Publications**

Basic Information

Observed MW

67kDa

Calculated MW

68kDa

Category

Polyclonal Antibodies

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Conjugation

Unconjugated

Background

Predicted to enable histone deacetylase binding activity; transcription coactivator activity; and transcription corepressor activity. Involved in negative regulation of transcription, DNA-templated. Located in nucleoplasm. Part of NuRD complex.

Recommended Dilutions

WB 1:500 - 1:2000

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Immunogen Information

Gene ID

57504

Swiss Prot

Q9BTC8

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

MTA3

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of lysates from A-549 cells, using MTA3 Rabbit pAb (A2328) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 30s.

