

PE/Cyanine7 anti Human CD223/LAG-3 mAb

Catalog NO.: A27492

Basic Information

Observed MW

Calculated MW

57kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

PE-Cy7. Ex:565nm. Em:778nm.

Clone number

AMC0821

Background

Lymphocyte-activation protein 3 belongs to Ig superfamily and contains 4 extracellular Ig-like domains. The LAG3 gene contains 8 exons. The sequence data, exon/intron organization, and chromosomal localization all indicate a close relationship of LAG3 to CD4.

Recommended Dilutions

FC 5 µl per 10⁶ cells in 100 µl volume

Immunogen Information

Gene ID

3902

Swiss Prot

P18627

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

CD223

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Mouse

Isotype

IgG1, kappa

Purification

Affinity purification

Storage

Store at 4°C. Avoid freeze.

Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.2.

