

STAT3 Rabbit mAb

Catalog NO.: A27918 **Recombinant**

Basic Information

Observed MW

79 kDa/86 kDa

Calculated MW

88 kDa

Category

Monoclonal Antibodies

Applications

WB, IHC-P, ChIP, ELISA

Cross-Reactivity

Human, Mouse, Rat, Monkey

Conjugation

Unconjugated

Clone number

ARC70808

Recommended Dilutions

WB 1:5000 - 1:20000

IHC-P 1:1000 - 1:4000

ChIP 2 µg antibody for 5 µg-15 µg of Chromatin

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your specific assay requirements.

Contact

Background

The protein encoded by this gene is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. This protein is activated through phosphorylation in response to various cytokines and growth factors including IFNs, EGF, IL5, IL6, HGF, LIF and BMP2. This protein mediates the expression of a variety of genes in response to cell stimuli, and thus plays a key role in many cellular processes such as cell growth and apoptosis. The small GTPase Rac1 has been shown to bind and regulate the activity of this protein. PIAS3 protein is a specific inhibitor of this protein. This gene also plays a role in regulating host response to viral and bacterial infections. Mutations in this gene are associated with infantile-onset multisystem autoimmune disease and hyper-immunoglobulin E syndrome.

Immunogen Information

Gene ID	Swiss Prot
6774	P40763

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

APRF; HIES; ADMIO; ADMIO1

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.09% Sodium azide,0.05% BSA,50% glycerol,pH7.3.

Validation Data



Western blot analysis of various lysates using STAT3 Rabbit mAb (A27918) at 1:18000 dilution incubated at room temperature for 1.5 hours.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.


Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Negative control (NC): PC-3

Exposure time: 45s.



Western blot analysis of various lysates using STAT3 Rabbit mAb (A27918) at 1:10000 dilution incubated at room temperature for 1.5 hours.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25 µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.


Detection: ECL Basic Kit (RM00020).

Exposure time: 1s.




Chromatin immunoprecipitation was performed with 15 μ g of cross-linked chromatin from HeLa cells were treated with IFN- α (100 ng/mL, 30 min), using 2 μ g of STAT3 Rabbit mAb (A27918) and Rabbit IgG isotype control (AC042). The enrichment of immunoprecipitated DNA at different genomic loci was examined by quantitative PCR. The histogram compares the ratio of the immunoprecipitated DNA to the input at given loci.

Validation Data



Immunohistochemistry analysis of paraffin-embedded Human colon carcinoma tissue using STAT3 Rabbit mAb (A27918) at a dilution of 1:1000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using STAT3 Rabbit mAb (A27918) at a dilution of 1:1000 (40x lens). High pressure antigen retrieval performed with 0.01M Tris-EDTA Buffer (pH 9.0) prior to IHC staining.

