

Biotin Rabbit anti-Mouse IgA (Fc) mAb

Catalog NO.: A28356

Basic Information

Observed MW

Calculated MW

Category

Monoclonal Antibodies

Applications

FC (intra)

Cross-Reactivity

Mouse

Conjugation

Biotin

Clone number

ARC71353

Background

Recommended Dilutions

FC (intra)

≤0.25 µg per million cells in 100 µl volume

Immunogen Information

Gene ID

Swiss Prot

/

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS with 0.09% Sodium azide,0.2% BSA,pH7.3.		

Validation Data



Flow cytometry: 1X10⁶ 293T (Transfection) cells were intracellularly-stained with Biotin Rabbit IgG isotype control (A25626,5 µl/Test,left) or Biotin Rabbit anti-Mouse IgA (Fc) mAb (A28356,0.25 µg,right), followed by APC Streptavidin staining.

