

Rabbit anti-Human CD7(scFv) mAb

Catalog NO.: A28936

Basic Information

Observed MW

Calculated MW

25 kDa

Category

Monoclonal Antibodies

Applications

FC

Cross-Reactivity

Human

Conjugation

Unconjugated

Clone number

ARC54668 (scFv)

Background

This gene encodes a transmembrane protein which is a member of the immunoglobulin superfamily. This protein is found on thymocytes and mature T cells. It plays an essential role in T-cell interactions and also in T-cell/B-cell interaction during early lymphoid development.

Recommended Dilutions

FC 1:100 - 1:500

Immunogen Information

Gene ID

924

Swiss Prot

P09564

Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

Synonyms

GP40; TP41; Tp40; LEU-9


Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source	Isotype	Purification
Rabbit	IgG	Affinity purification
Storage		
Store at 2-8°C. Avoid freeze.		
Buffer: PBS,pH7.3,2mM DTT		

Validation Data



Flow cytometry: 1X10⁶ Human PBMC were surface-stained with Rabbit IgG isotype control (AC042,2 µg/mL,left) or Rabbit anti-Human CD7(scFv) mAb (A28936,2 µg/mL,right), followed by Alexa Fluor® 647 His tag staining. Cells in the lymphocyte gate were used for analysis.

