

# SOS1 Rabbit pAb

Catalog NO.: A3272 **1 Publications**

## Basic Information

### Observed MW

175kDa/200kDa

### Calculated MW

152kDa

### Category

Polyclonal Antibodies

### Applications

WB, ELISA

### Cross-Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

## Background

This gene encodes a protein that is a guanine nucleotide exchange factor for RAS proteins, membrane proteins that bind guanine nucleotides and participate in signal transduction pathways. GTP binding activates and GTP hydrolysis inactivates RAS proteins. The product of this gene may regulate RAS proteins by facilitating the exchange of GTP for GDP. Mutations in this gene are associated with gingival fibromatosis 1 and Noonan syndrome type 4.

## Recommended Dilutions

**WB** 1:500 - 1:2000

### ELISA

Recommended starting concentration is 1 µg/mL.  
Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

| Gene ID | Swiss Prot |
|---------|------------|
| 6654    | Q07889     |

### Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

### Synonyms

GF1; HGF; NS4; GGF1; GINGF; SOS-1; SOS1


## Contact

## Product Information

|   |  |                      |
|---|--|----------------------|
| ☒ |  | order@abclonal.com   |
| ☒ |  | service@abclonal.com |
| ☒ |  | www.abclonal.com     |

|                                                                                                                                       |                |                       |
|---------------------------------------------------------------------------------------------------------------------------------------|----------------|-----------------------|
| <b>Source</b>                                                                                                                         | <b>Isotype</b> | <b>Purification</b>   |
| Rabbit                                                                                                                                | IgG            | Affinity purification |
| <b>Storage</b>                                                                                                                        |                |                       |
| Store at -20°C. Avoid freeze / thaw cycles.                                                                                           |                |                       |
| Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3. |                |                       |

Validation Data



Western blot analysis of various lysates using SOS1 Rabbit pAb (A3272) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.

