

Epac1/RAPGEF3 Rabbit mAb

Catalog NO.: A4149

Recombinant

1 Publications

Basic Information

Observed MW

100kDa

Calculated MW

104kDa

Category

Monoclonal Antibodies

Applications

WB, IHC-P, ELISA

Cross-Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Clone number

ARC0908

Background

Enables guanyl-nucleotide exchange factor activity and protein domain specific binding activity. Involved in several processes, including positive regulation of protein modification process

Recommended Dilutions

WB 1:1000 - 1:2000

IHC-P 1:100 - 1:400

ELISA

Recommended starting concentration is 1 µg/mL.
Please optimize the concentration based on your
specific assay requirements

Immunogen Information

Gene ID

10411

Swiss Prot

095398

Immunogen

Synthetic peptide. This information is considered to be commercially sensitive.

Synonyms

EPAC; Epac1/RAPGEF3

Contact

Product Information

☒		order@abclonal.com
☒		service@abclonal.com
☒		www.abclonal.com

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS containing 50% glycerol and 0.05% BSA, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3.

Validation Data



Western blot analysis of various lysates using Epac1/RAPGEF3 Rabbit mAb (A4149) at 1:1000 dilution.


Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 60s.



Immunohistochemistry analysis of paraffin-embedded Human placenta using Epac1/RAPGEF3 Rabbit mAb (A4149) at dilution of 1:100 (40x lens). Microwave antigen retrieval performed with 0.01M PBS Buffer (pH 7.2) prior to IHC staining.

