

# Langerin/CD207 Rabbit pAb

Catalog NO.: A6553

## Basic Information

### Observed MW

40-50kDa

### Calculated MW

37kDa

### Category

Polyclonal Antibodies

### Applications

WB, ELISA

### Cross-Reactivity

Human, Mouse, Rat

### Conjugation

Unconjugated

## Background

The protein encoded by this gene is expressed only in Langerhans cells which are immature dendritic cells of the epidermis and mucosa. It is localized in the Birbeck granules, organelles present in the cytoplasm of Langerhans cells and consisting of superimposed and zippered membranes. It is a C-type lectin with mannose binding specificity, and it has been proposed that mannose binding by this protein leads to internalization of antigen into Birbeck granules and providing access to a nonclassical antigen-processing pathway. Mutations in this gene result in Birbeck granules deficiency or loss of sugar binding activity.

## Recommended Dilutions

**WB** 1:500 - 1:1000

### ELISA

Recommended starting concentration is 1 µg/mL.  
Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

50489

### Swiss Prot

Q9UJ71

### Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

### Synonyms

CLEC4K; Langerin/CD207


## Contact

## Product Information

|                                     |  |                      |
|-------------------------------------|--|----------------------|
| <input checked="" type="checkbox"/> |  | order@abclonal.com   |
| <input checked="" type="checkbox"/> |  | service@abclonal.com |
| <input checked="" type="checkbox"/> |  | www.abclonal.com     |

|                                                                                                                                       |                |                       |
|---------------------------------------------------------------------------------------------------------------------------------------|----------------|-----------------------|
| <b>Source</b>                                                                                                                         | <b>Isotype</b> | <b>Purification</b>   |
| Rabbit                                                                                                                                | IgG            | Affinity purification |
| <b>Storage</b>                                                                                                                        |                |                       |
| Store at -20°C. Avoid freeze / thaw cycles.                                                                                           |                |                       |
| Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3. |                |                       |

Validation Data



Western blot analysis of various lysates using Langerin/CD207 Rabbit pAb (A6553) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.

