

# SAE1 Rabbit pAb

Catalog NO.: A9960

## Basic Information

### Observed MW

38kDa

### Calculated MW

38kDa

### Category

Polyclonal Antibodies

### Applications

WB, ELISA

### Cross-Reactivity

Human, Mouse

### Conjugation

Unconjugated

## Background

Posttranslational modification of proteins by the addition of the small protein SUMO (see SUMO1; MIM 601912), or sumoylation, regulates protein structure and intracellular localization. SAE1 and UBA2 (MIM 613295) form a heterodimer that functions as a SUMO-activating enzyme for the sumoylation of proteins (Okuma et al., 1999 [PubMed 9920803]).

## Recommended Dilutions

**WB** 1:500 - 1:2000

### ELISA

Recommended starting concentration is 1 µg/mL.  
Please optimize the concentration based on your specific assay requirements.

## Immunogen Information

### Gene ID

10055

### Swiss Prot

Q9UBE0

### Immunogen

Recombinant protein (or fragment).This information is considered to be commercially sensitive.

### Synonyms

AOS1; SUA1; UBLE1A; HSPC140; SAE1


## Contact

## Product Information

|                                     |  |                      |
|-------------------------------------|--|----------------------|
| <input checked="" type="checkbox"/> |  | order@abclonal.com   |
| <input checked="" type="checkbox"/> |  | service@abclonal.com |
| <input checked="" type="checkbox"/> |  | www.abclonal.com     |

|                                                                                                                                       |                |                       |
|---------------------------------------------------------------------------------------------------------------------------------------|----------------|-----------------------|
| <b>Source</b>                                                                                                                         | <b>Isotype</b> | <b>Purification</b>   |
| Rabbit                                                                                                                                | IgG            | Affinity purification |
| <b>Storage</b>                                                                                                                        |                |                       |
| Store at -20°C. Avoid freeze / thaw cycles.                                                                                           |                |                       |
| Buffer: PBS containing 50% glycerol, preserved with proclin300 or sodium azide (as specified on the Certificate of Analysis), pH 7.3. |                |                       |

Validation Data



Western blot analysis of various lysates using SAE1 Rabbit pAb (A9960) at 1:1000 dilution.

Secondary antibody: HRP-conjugated Goat anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.

Lysates/proteins: 25µg per lane.

Blocking buffer: 3% nonfat dry milk in TBST.

Detection: ECL Basic Kit (RM00020).

Exposure time: 90s.

