

Human/Monkey Treg Panel

Catalog Number: RK40022

Size: 1 Box

◆ Description

The 7-fluorescence Human/Monkey Treg Panel has been designed for identify and phenotype T regulatory cells in Human/Monkey blood or tissue samples.

ABflo® multi-fluorescence panel RK40022 covers the core CD markers of CD45, CD3, CD4, CD8a, CD14, CD25 and FOXP3. Our plots showed clear subgroups of CD14⁺CD3⁺ T lymphocytes, CD3⁺CD4⁺ T lymphocytes, CD3⁺CD8⁺ T lymphocytes and CD4⁺CD25⁺FOXP3⁺ Regulatory T cells (cite figures).

◆ Applications

FC

◆ Reactivity

Human/Monkey

◆ Components

Target	Catalog No	Fluorochrome	Ex/Em (nm)	Isotype	Clone	Size	Vol. Per Test
CD45	A27932	ABflo® 700	690/713	Rabbit IgG	ARC5024	100T	5 µL
CD3	A27111	APC	648/664	Rabbit IgG	ARC68274	100T	5 µL
CD4	A26523	ABflo® 450	406/445	Rabbit IgG	ARC5142	100T	5 µL

CD8a	A26731	ABflo [®] 500	410/501	Rabbit IgG	ARC55248	100T	5 µL
CD14	A25071	ABflo [®] 488	491/516	Rabbit IgG	ARC65657	100T	5 µL
CD25	A27400	PE/Cyanine7	565/778	Rabbit IgG	ARC58045	100T	5 µL
FOXP3	A27686	PE	565/574	Rabbit IgG	ARC74142	100T	5 µL


Note:

This reagent has been pre-diluted for use at the recommended Volume per Test. We typically use 1×10^6 cells in a 100 µl experimental sample (a test).

◆ Storage

Store at 4°C. Avoid freeze/thaw cycles.

◆ Validation Data



Multicolor flow cytometry panel to identify Human/Monkey Treg cell subsets. Cynomolgus monkey blood were surface stained with ABflo[®] 700 Rabbit anti-Human/Monkey CD45 mAb (A27932, 5 µL/Test), APC Rabbit anti-Human/Monkey CD3 mAb (A27111, 5 µL/Test), ABflo[®] 450 Rabbit anti-Human/Monkey CD4 mAb (A26523, 5 µL/Test), ABflo[®] 500 Rabbit anti-Human/Monkey CD8a mAb (A26731, 5 µL/Test), ABflo[®] 488 Rabbit

anti-Human/Monkey CD14 mAb (A25071,5 µL/Test), PE/Cyanine7 Rabbit anti-Human/Monkey CD25 mAb (A27400,5 µL/Test), PE Rabbit anti-Human/Mouse/Monkey FOXP3 mAb (A27686,5 µL/Test). All antibodies are validated for Flow cytometry. CD3⁺ T lymphocytes are defined as CD14⁻CD3⁺; CD4⁺ T lymphocytes are defined as CD3⁺CD4⁺; CD8⁺ T lymphocytes are defined as CD3⁺CD8⁺; Regulatory T cells are defined as CD4⁺CD25⁺FOXP3⁺.